

plastic solutions

FOR MEDICAL DEVICES



Seraclean® (Serumfilter for blood analysis)

GEMÜ in corporation with:

KUNSTSTOFF TRENDS
MÄRKTE · PRODUKTE · ANWENDUNGEN

BIOforum
FORSCHUNG · FACHKOLLEGEN · SERVICE

**STERIL
TECHNIK**

Management &
Krankenhäuser

GIT VERLAG
A Wiley Company



GEMÜ® MEDICAL PLASTICS
PROCESSING DIVISION

REPRINT

Cost Efficacy and Value of Time-tested Systems in Clinical Diagnostics

Lars Rominger

Several times in the recent past, heads of central laboratories and those in charge of purchasing and procurement in clinics have been faced with the decision as to whether they should invest in new medical technologies. In most instances the criteria for decision-making are clear, demonstrable, and can be evidenced in terms of economic efficiency as well as the medical value of the technology. However, in some instances technological advances are the only driving force for those offering new products in this sector. For laboratory executives this means that a clear system-based margin must be drawn between new investments in uncertain products on the one hand, and the stable, safe and continued application of time-tested systems and processes on the other.



About the author:

Lars Rominger is Head of Marketing/Sales, Medical Engineering, at GEMÜ GmbH Switzerland.

For several years now there has been a quest for a time-tested, safe, medically stable, financially established and economical principle in medical diagnostics. Human blood for analysis is centrifuged. By this process blood is fractionated into two pha-

ses: the blood clot remains below and the serum is collected above. After centrifugation a blood filter is introduced into the centrifuge glass and the serum is filtered. This standardized procedure is used to filter out any residual clots or other components that may still be present in the fluid above.

The purpose of this filtering procedure is to avoid any risk of falsifying the test result and also prevent sedimentation in the instruments used for analysis.

Gel versus Blood Filter

This standardized procedure is being applied successfully for several years and is now established as a standard in nearly all laboratories. The system is time-tested and absolutely reliable. The use of blood filters produces reliable test substances for the subsequent measurement. The company that has been manufacturing this product so far recently informed its market partners that the above mentioned blood filter will no longer be pro-

duced or distributed. It was also stated that the manufacturing companies, after several years of intensive and laborious research, are now backing a procedure that will gradually replace the previous system and also simplify it. A principle based on a gel procedure was developed and the product manufactured according to this principle. Separation gel tubes are in use for a few years now. They were developed to improve the process of pre-analysis. The production of the blood filters was stopped because the manufacturing company presumed that the increasing use of separation gel tubes would replace the blood filters. The manufacturer is actively working on marketing and introducing the new gel-based procedure. However, according to users in central laboratories, the market maturity of this procedure is highly disputed, although the technology of analysis instruments has made significant progress.

As the production and delivery of the manufacturer's time-tested and esteemed blood filters has been suspended, laboratory executives were now



Seventy staff members are employed by GEMÜ GmbH (private limited company) in Switzerland - a competence center for plastics and medical technology, part of the GEMÜ group which operates on a worldwide basis, with its head office in the town of Ingelfingen-Criesbach close to Stuttgart (Germany). The GEMÜ enterprise is active in the fields of medicine, medical technology, the pharmaceutical and food industry, process engineering and microchip production. Processing technology at GEMÜ GmbH in Switzerland focuses on the injection molding of plastics. In this sector the company offers system-based solutions for plastics in the medical industry, in biotechnology, the semiconductor and electronics industry. At the Swiss office in Rotkreuz, extensive resources have been invested over the last few years in clean room concepts for the production of medical-technical products. GEMÜ engages more than 850 staff members throughout the world.

Re-launching of Time-tested Products

Based on a market analysis, the managers of GEMÜ GmbH in Switzerland recognized the precarious situation early and developed measures to resolve this critical state of affairs. Within a short period of time, a task force under Lars Rominger was established within the company and empowered with the authorities required to do the job. In a few weeks the task force designed an innovative system-based solution to the problem, based on

stringent process-oriented project management. This was the driving force behind the newly developed serum filter **Seraclean** - a solution to fulfill a specific need, offered at a price in line with market conditions.

The product specifications provide evidence of GEMÜ GmbH's (Switzerland) extensive specialized know-how in the field of system-based solutions in plastics for medicine and the industry.

The Products and how they Work

Human blood for analysis is centrifuged. By this procedure blood is fractionated into two phases: the blood clot remains below and the serum is collected above. After centrifugation, a blood filter is inserted into the centrifuge glass and the serum filtered by this procedure. The purpose of this standardized process is to filter out any residual clots or other components in the fluid above.

Integrated Additional Technical Value

In contrast to other filter systems Seraclean is equipped with a sealing-lip that never ceases to function because of the specific tool manufacturing technology and the manufacturing processes used to produce this instrument. Besides, based on a number of detailed technical discussions with laboratory executives and users of serum filters, Seraclean has been designed such that it provides maximum functionality while its application is as user-friendly as possible.



being forced to focus their attention and interest exclusively on the gel-based procedure. In the current environment of cost-sensitivity and economic pressures in the public health sector, requests for investments in new devices and/or systems are received with little enthusiasm by the financial departments of medical clinics.

The situation is rendered more difficult by the fact that personnel - in this case medical laboratory staff - have to be re-trained within the hospital. In addition to new investments, this change involves additional effort in terms of training time and monetary resources for re-training materials.

In every case, decision-makers among the medical staff are mainly concerned with maintenance of the service unit, work processes, and the required purity and reliability of blood samples. One of the aspects of the core competence of a laboratory is its ability to deliver a large number of absolutely reliable blood analysis within a short period of time. This requirement must be fulfilled. It poses a challenge, particularly for medical technology.



Market Acceptance from the 0 Series Onward

In numerous conversations and detailed hearings with clinical laboratory executives, both the task force and the executives of GEMÜ received adequate confirmation of the necessity and appropriateness of this developmental project. A re-launch of a long-standing time-tested serum filter on a high-tech synthetic basis, in line with the highest standards of clean room technology, is a high-level market opportunity. Market partners at all levels of the hierarchy and in all positions have asked for a simple and safe product in alignment with real market conditions. The acceptance of the newly introduced **Seraclean** filter in the market is an undisputed fact. This was confirmed in several conversations that Lars Rominger held with Dr. Engler and Ms. Bossart from Kantonsspital St. Gallen (Institute for Clinical Chemistry and Hematology), and Lukas Bestmann, Head of the Department of General Analytics from the University Hospital of Zurich (Institute for Clinical Chemistry).

The Direct Advantage for the Client is Self-evident:

- The instrument requires no new financial investments in central laboratories. It involves no costly revisions of specialized devices due to gel-induced sediments or obstructions. Seraclean works according to a well established and well known work process.
- In order to achieve reliable test results with the GEL procedure as well, Seraclean may be used as an additional process.
- Seraclean has none of the disadvantages of the GEL procedure: disturbances in the instrumental analysis due to infinitesimal quantities of GEL entering the serum; undesired post-coagulation; the gel is not inert and therefore sensitive to cold.
- Gels absorb large molecules (e.g. medications) which may falsify the test results.
- In contrast to other serum filter procedures Seraclean is very user-friendly and fulfills the highest standards of functionality.

Expert

Interview with the specialists Dr. Hanna

Engler from Kantonsspital St. Gallen

(Institute for Clinical Chemistry and He-

matology) and Lukas Bestmann, Head

of the Department of General Analytics

from the University Hospital of Zurich

(Institute for Clinical Chemistry).

Can you say a few words about your experience with the gel-based procedure?

Dr. Engler: Separation gels simplify the work process because they create a division - a layer - between the blood clot and the serum/plasma during the centrifugation step. This layer prevents interaction between the cell components and the blood clot. Besides, the serum/plasma layer can be transferred into a tube (which is then used to transport the sample) by means of simple decantation. This is particularly advantageous for samples which have to be sent to the clinical-chemical laboratory by mail.

The disadvantages of separation gel tubes include the fact that some serum samples may post-coagulate after the centrifugation step, and also gel particles may be released occasionally. If the phenomenon remains undetected it may lead to false test results. In extreme cases these gel particles may cause obstructions in the autoanalyzer system, involving significant costs later on. It should also be remembered that separation gels may not be inert. Thus, the analytes contained in the sample

Opinions

may be re-absorbed. This could lead to falsely low test results and inappropriate clinical decisions.

Mr. Bestmann: Gels absorb many large molecules such as medications. The values yielded by the analysis may be 40% below the actual value. This does not happen with a filter system. Such false low values can be dangerous. They may distort the general impression – they would falsely indicate that the patient has very low values of a medication in his/her blood; the doctor may increase the dose and the patient may reach a critically high value (in the worst case this could be in the toxic range). Thus, it will be very difficult to steer pharmacological therapy. Besides, gels create technical problems in special devices such as GC-MS, LC-MS, MALDI-TOF etc. Gels may also interfere with certain ELISA's, clog devices, and lead to costly revisions.

What, in your opinion, are the advantages of the serum filter procedure?

Dr. Engler: Serum filters serve the same purpose as separation gels. They introduce a barrier between the blood clot and the serum/plasma, which prevents pre-analytical problems due to interactions between cell components and serum/plasma. Thus the original tubes can be used in the analytical devices; this permits positive identification of the sample.

As these serum filters are introduced in the sampling tube after the centrifugation step, it is ensured that post-coagulations and other floating particles are removed from the serum/plasma sample. As serum filters can be used with all types of blood sampling tubes, the pipetting step is not required if the blood sampling tube is used without a separating aid.

Mr. Bestmann: Subsequent orders can be made for a longer period of time because the analyses are falsified more slowly. One example is LDH which, in the absence of a filter, is diffused from red blood cells into plasma after a few hours and

generates false high values. When a filter is used, the enzyme is withheld and the clinic may send a subsequent order within 12 hours; this could spare the patient the discomfort of additional blood sampling.

What, in your opinion, are the specific advantages of Seraclean?

Dr. Engler: Seraclean is introduced in the sample container by means of pressure. The long tube allows the filter to be used for every insertion depth. The material is not brittle. Thus, protruding parts of the filter can be cut away so that the original containers of samples with a very high hematocrit content (large percentage of blood cells and relatively low quantity of serum/plasma) can be used directly in the analyzer even after the filter has been introduced. Another advantage is the rounded upper rim and the hole in the upper end. These reduce the risk of injury and the filter can be pushed directly into the container with the thumb because the air is able to escape during the insertion procedure.

Mr. Bestmann: Seraclean has no ascending pipe – this makes the filter suitable for the analysis instruments. In devices with an ascending pipe, inexact adjustment may damage pipette needles, and pipette needles are very expensive. As Seraclean is similar to its predecessor, it is easy to handle. Thanks to the long neck of the filter, it can be inserted to the appropriate depth (to the margin of the blood clot) for each sample. The hole at the upper end allows the filter to be pushed downward with the thumb; this causes air to escape during the insertion procedure. Some devices lack this advantage.

What are your expectations of a product that is part of the instrumentation used for clinical diagnostic procedures?

Dr. Engler: Simplicity, reliability, efficiency, no interaction with analytes or measuring systems, and low extra costs. It would be desirable if the

serum filter could be integrated into pre-analytical systems. For instance, automatic insertion and adjustment of the height of the filter in the blood sampling tube. This should occur after the automated centrifugation step.

Mr. Bestmann: Reliability, suitability, simplicity, efficiency – preferably, there should be an automatic system of pressing the filter into the blood withdrawal system.

Dr. Engler and Mr. Bestmann – thank you very much for this interesting and informative interview.



Dr. Hanna Engler Lukas Bestmann

Contact:

GEMÜ GmbH (Switzerland)
Plastics and Medical Technology
Lars Rominger
Head of Marketing/Sales, Medical Engineering
P.O. Box
Lettenstrasse 3
CH-6343 Rotkreuz
Tel. +41 (0)41 799 05 05
Fax +41 (0)41 799 05 85
Lars.Rominger@gemue.ch
www.gemue.ch



Here you see the strongest argument
for holding on to proven serum filters:
Seraclean®



Why do so many leading pharmacists opt for Seraclean®? Perhaps because we have taken up a proven, unmatched serum-filtering process. Probably because we have further perfected it in cooperation with renowned laboratories. Most certainly, however, because they receive a product which, in terms of safety, functionality and user-friendliness, sets new standards and, in addition, spares them expensive new investments. Seraclean® – and you will have the best arguments in hand.



Purity your requirement our way of living



Seraclean®



System solutions in plastics for medicine, pharmacy and industry. GEMÜ Switzerland provides complete system solutions in plastics for the medical device industry and industrial applications ranging from engineering and mold manufacturing, project and quality management, rapid prototyping, mold-flow analysis to injection molding, manufacturing, cleaning, testing, final assembly and packaging, depending on customer requirements – Process-oriented from concept to finished product!



Have we aroused your interest?

Achieve Your Target With Us!

By Fax to GEMÜ +41 (0) 41 799 05 85

- Yes, we are interested in **plastics processing under clean-room conditions at GEMÜ**
- Yes, we are interested in **Seraclean® (Serumfilter for blood analysis)**

Please call us on Telephone No.: _____

Our company is active in the following field(s):

- | | |
|---|--|
| <input type="checkbox"/> Medical engineering | <input type="checkbox"/> Automobile industry |
| <input type="checkbox"/> Pharmaceuticals industry | <input type="checkbox"/> LED production |
| <input type="checkbox"/> Biotechnology | <input type="checkbox"/> Solar-cell manufacture |
| <input type="checkbox"/> Foodstuffs industry | <input type="checkbox"/> Aviation and space technology |
| <input type="checkbox"/> Chemical production | <input type="checkbox"/> Telecommunications |
| <input type="checkbox"/> Semi-conductor industry | <input type="checkbox"/> others |

This Fax is from:

Address: _____

Date/Signature: _____

Your contact at GEMÜ: